

Checking The Defibrillator

1 Does it have a cabinet

Various questions about the cabinet. All OK.



2 Is it ready to be used

Press the middle button (red box in picture)

Shows battery and pad dates



3 Is the battery inserted and OK



4 Are the main pads connected and in date



5 Is the rescue kit ready for use

Check pack is behind main defibrillator



6 Complete report on website



<https://www.defibcheck.org/>

Website to file report

[Simon - simon.kearsley@symmetry.co.uk](mailto:simon.kearsley@symmetry.co.uk)

Scheme Coordinator

Michael - clerk@hankerton-pc.gov.uk

Parish Council Clerk

Jaqui Erskine-Crum - dolmansfarm@aol.com

Parish Council Chair